

Command:RESET (~)

Shortcut	~
Format	RESET

Additional Notes about usage on Tiger

RESET is usually a " Broadcast command" but can be used as a Card-Addressed Command as well. When addressed to a specific card, it resets that card and then the Comm card (required to communicate with the reset card) without affecting other installed cards.

<html>
</html> This command causes the controller to do a software reset. A software reset reinitializes all variables back to their pre-assigned values.

A reset is different from a power cycle. This command does not save positions, commands like [SETHOME](#) require a power cycle instead of a reset.

Example

```
RESET  
:A
```

[commands](#), [tiger](#), [ms2000](#)

From:

<https://docs.asiimaging.com/> - **Applied Scientific Instrumentation**

Permanent link:

<https://docs.asiimaging.com/commands/reset?rev=1603927929>

Last update: **2021/09/23 15:23**

