

MATLAB

We currently don't have any official support for MATLAB.

You can still use MATLAB by sending serial commands to our controllers. See our [Quick Start on Serial Commands](#) page.

Here is a unofficial script to control ASI stages using an MS2000:

<https://www.mathworks.com/matlabcentral/fileexchange/19428-asi-microscope-stage-mover>.

You can utilize Micro-Manager through Matlab. Micro-Manager provides a uniform API to control many microscopy devices and ASI's controllers are very well supported. See

https://micro-manager.org/wiki/Matlab_Configuration.

[matlab](#)

From:

<https://docs.asiimaging.com/> - **Applied Scientific Instrumentation**

Permanent link:

<https://docs.asiimaging.com/matlab?rev=1632431711>

Last update: **2021/09/23 17:15**

