

Microinjector Pipette Holder MIMPH-MPIP Kit

The Microinjector Pipette Holder is designed to hold glass drawn micropipettes in a stable, easy to mount, easy to clean, autoclaveable, stainless steel body.

Warning! Micro drawn glass pipettes pose a cut or puncture risk. Always secure the holder away from hand motion and do not let the pipette point at yourself or another person. Always secure the glass capillary pipette carefully in the collet; an unsecured pipette can shoot out of the assembly at high velocity due to high pressures from the injection device.

[Microinjector Pipette Holder MIMPH-MPIP Kit](#) - PDF Manual

[mppi](#), [manual](#)

From:

<https://docs.asiimaging.com/> - **Applied Scientific Instrumentation**

Permanent link:

https://docs.asiimaging.com/mimph_mpip_kit?rev=1732664606

Last update: **2024/11/26 18:43**

